Product Description

Vergenix™STR is an advanced soft tissue repair device and a biodegradable scaffold made of lyophilized Type I rhCollagen and calcium chloride. When the device is combined with autologous Platelet Rich Plasma (PRP), it is injected and forms a clot in the vicinity of the tendon injury site. Vergenix™STR is intended for the treatment of tendinopathy (tendon inflammation). Vergenix™STR serves as a scaffold for fibroblast recruitment. Vergenix™STR contains collagen that forms a 3 ml suspension following mixing with PRP. The device is supplied sterile in a single use kit, packaged in Tyvek sealed blister containing one syringe with rhCollagen and calcium chloride, one empty syringe, one luer lock connector, and one 19G needle.

Indications for Use

Vergenix™STR is indicated for the treatment of tendinopathy.

Precautions

The device is for single use only. Vergenix™STR injection together with PRP is intended for use by trained orthopedic surgeons in hospitals, clinics or private practices. Prior to the treatment with Vergenix™STR, a local anesthesia will be injected into the injection site. The blood collection for PRP separation, prior to mixing with Vergenix™STR, must be performed by a health care professional, such as a nurse or physician following training on the use of the commercially available PRP kit. Vergenix™STR must be used with autologous PRP only. That is, PRP must be injected into the same patient from whom the blood sample for PRP preparation was drawn to prevent contamination.

The PRP fraction can be used up to time limits indicated by the PRP kit collection product instructions for use, after drawing blood from patient. Vergenix™STR will be used following tendinopathy diagnosis by an orthopedic surgeon.

All blood collection aids and injection needle should be handled and disposed of with caution according to established policies of the health care setting. All injections and blood collection steps should be performed on a sanitized site. Do not use Vergenix™STR if foreign matter is present. Do not use Vergenix™STR if packaging is damaged. Do not re-sterilize Vergenix™STR and do not attempt to clean the product. Do not use a needle size that is smaller than 19G.

Note that the degree of accuracy of the 5 ml syringe is +/-0.2 ml.

Vergenix™STR must not be reused. Reuse of product will cause loss of sterility and risk of infection.

Contraindications

Vergenix™STR must not be used in patients with known sensitivity to collagen. Discontinue use if signs of allergy/hypersensitivity appear.

Adverse Reactions

The following events may potentially occur as with any type of injection: potential damage to nerve or blood vessels, hematoma, pain, edema formation and superficial infection.

Directions for Use

Note: Failure to follow Direction for Use could result in delayed, poor or no tendon healing.

Step 1 - Local anesthesia injection into the injured tendon

- Sanitize injection site.
- Inject 1ml Lidocaine or Marcaine to the site, using 25G needle.

Step 2 - PRP collection

- Follow the commercially available PRP kit instructions for use in order to isolate the PRP fraction.
- Keep PRP aside until used as per PRP kit instructions for use.

Step 3 - Vergenix™STR mixing with PRP - This action should be performed shortly before injection into the tendon

- Withdraw 2.7 ml of PRP to an empty 5cc syringe.
- Connect the luer-lock connector to the syringe containing PRP.
- Advance the syringe plunger to fill the luer-lock connector with PRP.
- Connect the syringe containing the dry collagen to the other side of the luer-lock.

Make sure that the injection area is ready, including sanitation and anesthesia. The collagen needs to be mixed with PRP shortly before the injection step.

- Transfer the PRP into the syringe containing dry collagen. Tap the loaded syringe on the palm of the hand, by holding the empty connected syringe ~20 times.
- Disconnect the empty syringe and remove air bubbles by injecting the PRP-collagen suspension until filling the Luer-lock connector.
- Reconnect the empty syringe and mix the solution by transferring it from one syringe to the other ~20 times.
- Tap again the loaded syringe on the palm of the hand, by holding the empty connected syringe ~20 times.
- Mix the material again between the two syringes.
- The suspension is ready when it has a uniform and a pale orange color.
- Transfer the suspension to one syringe.
- Disconnect the luer-lock connector and empty syringe.
- Open carefully the appropriate sterile needle and connect it to the syringe.

Step 4 - Vergenix™STR injection into the injured tendon

- Inject the suspension to the injured tendon.

Storage

Vergenix™STR should be stored at room temperature (5°C - 30°C). Do not use after expiration date. Device is designated for a single use only.

Do not use if individual pack is opened or appears damaged.

For Single Use Only.

DO NOT RESTERILIZE.

Do not use after expiration date.

Vergenix™STR is supplied sterile. Method of sterilization is Ethylene Oxide.

Date of issue: April 2019

Part Number: LB11 Rev 5

Manufacturer:

ColiPlant, Ltd.
4 Oppenheimer St., P.O.B 4132
Rehovot 7670104, Israel

EC Authorized Representative

Donawa Lifescience Consulting Srl
Piazza Albania, 10 | 00153 Rome, Italy

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CollPlant, Ltd.